

## LESS PAIN. FAST RECOVERY.

Technicians can precisely target a surface area by adjusting the energy levels delivered into the tissue from 0.5mm to 3.5mm and prompting dermal coagulation using the advanced penetration mechanism.

## OPTIMAL & UNIFORM THERMOLYSIS

An optimal amount of space between the microneedles to uniformly distribute RF energy throughout the target area and evenly coagulate without the risk of interference caused by nearby, overlapping RF energy.

## PRECISE & SAFE TREATMENT

Designed with an Enhanced Delay function that helps eliminate post-inflammatory hyperpigmentation (PIH) by safely delivering the precise amount of heat necessary to stimulate the growth of collagen in the tissue.

## FAST & CONVENIENT OPERATION

Easy to operate and user-friendly, with a 15-inch touch interface in full color. The system has various inspection and cleaning modes. Plus, the guided-user interface provides an avatar system to apply a precise energy level and penetration depth to the patient's tissue for proven results.

# NOUVELLE

#### MICRONEEDLE RF


The NOUVELLE Microneedle RF system facilitates non-surgical skin tightening using bi-polar RF energy and minimally invasive microneedles to drive heat deep within the patient's skin layers. The heat stimulates and remodels collagen and prompts the elastin cells to regenerate. The treatment requires little to no downtime and results in a skin surface that is firmer, smoother, and more even.









### **HOW IT WORKS**




Microneedle cartridge makes perfect contact to the skin.



Microneedles smoothly penetrate the skin with less pain.



Bipolar RF energy denaturalizes the tissue surrounding the microneedles.



Collagen regeneration and new elastin production process begins.

#### **TONE & TEXTURE**







**AFTER** 

#### **TEXTURE IMPROVEMENT**



**BEFORE** 



**AFTER** 

#### **SKIN TIGHTENING**



**BEFORE** 



**AFTER** 

#### **NECK LIFTING**



**BEFORE** 



**AFTER** 

#### **ACNE SCAR IMPROVEMENT**



**BEFORE** 



**AFTER** 

### **INDICATIONS**

Treats a variety of skin conditions and skin types, with very low risk of burns and post inflammatory hyperpigmentation, to combat visible signs of aging on the face and body.

- Skin Tightening
- Wrinkle Reduction
- Acne Scar Improvement
- Skin Tone & Texture Improvement
- Pore Reduction
- Scar Reduction
- Stretch Mark Reduction
- Double Chin Reduction
- Jowl Lifting
- Eyebrow Lifting
- Neck Lifting



#### **NOUVELLE MICRONEEDLE RF System Specifications**

Power	Frequency	Level	Treatment Time	Interval	User Interface	Cooling
Max 25W	2 MHz (Bipolar)	Levels 1-10 *Power(W)=VRmzx- IRms2xxR=VRms2/R	50 ms - 200 ms (10 ms intervals)	Single 1.0 second 1.5 second 2.0 second	15" Touch Screen	Air

#### **Technical Specifications**

Electrical Rating: Single Phase 110-230VAC; 60/60 Hz; Power Consumption 500 VA (Fuse: 250V/6.3A)

Dimensions Including Handpiece Cable Hanger: 24 inches x 17 inches x 63 inches (with handpiece hanger)

Weight: 8 lb - 3.7 kg

### **ABOUT NYLO AESTHETICS**

NYLO Aesthetics is a leading distributor of aesthetic medical devices, with a portfolio of products offering best-in-class clinical results and patient experiences. NYLO globally sources innovative, state-of-the-art technology for U.S. medical offices, including dermatologists, plastic surgeons, and cosmetic practitioners, as well as premier medi-spas, and luxury hotel and day spas. With a large showroom and training facility in New York City, and offices in key markets nationwide, NYLO offers an immersive, hands-on experience.



65 West 36th Street | 10th FL | New York, NY 10018 800-887-1966 | contactus@nyloaesthetics.com www.nyloaesthetics.com